wellion® AMADEUS

Wellion AMADEUS Blutzuckerteststreifen

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EN Wellion AMADEUS Blood glucose test strips

IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLION AMADEUS BLOOD GLUCOSE METER USER'S GUIDE BEFORE USING WELLION GALILEO BLOOD GLUCOSE TEST STRIPS.

FOR QUESTIONS AND ASSISTANCE, PLEASE CONTACT THE AUTHORIZED REPRESENTATIVE. (FOR CONTACT DETAILS PLEASE REFER TO THE INFORMATION PROVIDED AT THE END OF THIS INSERT)

Intended Use

The Wellion AMADEUS Blood Glucose Test Strip is intended to be used with the Wellion AMADEUS Blood Glucose Meter to quantitatively measure glucose in venous whole blood, neonatal whole blood or capillary whole blood taken from fingertips, palm, or forearm. The Wellion AMADEUS Blood Glucose Test Strip is plasma-calibrated for easy comparison to lab results. It is intended for self-testing by people with diabetes and for professional use. It is not to be used for the diagnosis or screening for diabetes.

Introduction

The Wellion AMADEUS Blood Glucose Test Strip uses an enzyme to measure blood glucose. When blood touches the tip of the test strip, it flows into the reaction zone. The enzyme reacts with glucose in blood and produces electrical current. The meter measures the current. The meter

shows the test result within 8 seconds.

Reagent Composition

Each cm² of test strip contains:
[FAD] Glucose Dehydrogenase 13,14%
Electron Shuttle 38,16%
Non-reactive ingredients 48,7%

Warnings and Precautions

The Wellion AMADEUS Blood Glucose Test Strip is for use outside the body (IN VITRO diagnostic use).

[1] Do not use test strips after their expiration date. Please check the expiration date on the test strip vial. [2] Do not use strips that are wet or damaged. [3] Do not reuse strips. [4] If your test result is below 70 mg/dL (3,9 mmol/L) or above 250 mg/dL (13,9 mmol/L), perform a control solution test to make sure your system is working properly. Then repeat testing using fingertip blood. If the result is still very high or very low, contact your healthcare professional immediately. [5] If you have symptoms that are inconsistent with your test results and you have eliminated common errors as described in the user's guide, contact your healthcare professional immediately. [6] Never make major changes in your diabetes treatment program or ignore symptoms without consulting your physician. [7] Operation temperature for meter and test strip is between 5°C–45°C (41°F–113°F). [8] Incorrect results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic-hypersmolar state, with or without ketosis. [9] Do not test critically ill patients with blood glucose meters. [10] Incorrect results may occur in individuals who are dehydrated.

Storage and Handling

[1] Always close the vial cap tightly after removing a test strip. This avoids improper exposure to moisture and direct sunlight. [2] Store the test strips between 4°C – 30°C (39°F – 86°F) and 10% - 85% relative humidity. Do not freeze. [3] Store out of direct sunlight. [4] Unopened test strips are stable until the expiration date printed on the vial under specified storage conditions. [5] Use within 6 months after first opening. [6] Do not handle the test strip with wet or dirty hands. [7] Used strips are biohazardous materials and can transmit bloodborne disease. Dispose them according to local government regulations to avoid injury or contamination.

Sample Collection and Preparation

Testing must be performed immediately after the sample is obtained.

Test Procedure

Please refer to the chapter „Testing Your Blood Glucose“ in the Wellion AMADEUS Blood Glucose Meter User's Guide.

Alternate Site Testing

It is also possible to draw blood samples from palm and forearm. This is called alternate site testing (AST). AST results may differ from fingertip readings.¹⁰

AST should only be used:

- [1] In a pre-meal or fasting state (more than 2 hours last meal)
- [2] 2 hours or more after insulin application
- [3] 2 hours or more after exercise

Alternate Site Testing SHOULD NOT be used when:

- [1] You have hypoglycemic unawareness (not able to tell if you have low blood sugar).
- [2] Within 2 hours of a meal, exercise, or medication.
- [3] You operate machinery or drive a car.
- [4] You are sick.
- [5] You think your blood glucose is low.
- [6] Your AST results do not reflect the way you feel.
- [7] You suspect hyperglycemia.
- [8] Your blood glucose levels fluctuate frequently.
- [9] You are pregnant.

Consult your healthcare professional to decide if alternate site testing is right for you.

Expected Values

Consult with your physician or healthcare professional to determine an appropriate blood glucose target range for you.

Quality Control

Three different levels (Level 1, Level 2 or Level 3) of control solutions are available for quality checks. Follow the User's Guide instructions. It is recommended to perform control solution testing:

- [1] If your test results do not agree with how you feel.
- [2] Once per week to make sure that both the meter and the test strips are working properly.
- [3] If your test strips were stored at temperature and humidity outside specified storage conditions.
- [4] When you use your meter for the first time.
- [5] When you open a new vial of test strips.
- [6] To practice your testing technique.
- [7] If you dropped your meter.

Only use Wellion AMADEUS Control Solution. Your test results should be within the specified range provided on the test strip vial. Repeat control solution testing if results are outside specified ranges. Results outside the specified range may be due to:

- [1] Operating errors.
- [2] Expired or contaminated control solution
- [3] Damaged test strip
- [4] Meter malfunction

If the results repeatedly are outside the predefined range, contact the authorized representative whose contact details are provided at the end of this insert.

Troubleshooting

- [1] Confirm that the test strip is not expired.
- [2] Ensure proper blood sample application. If the blood sample volume is not sufficient, "Er3" is displayed. DO NOT add another drop of blood. Remove the test strip from the meter and retest with a new test strip.
- [3] Check the performance of the system with a control solution test.
- [4] Refer to the chapter "Troubleshooting" in the User's Guide for further support.

Limitations

- [1] DO NOT use plasma or serum samples.
- [2] DO NOT test on arterial blood samples.
- [3] Altitudes up to 3150m (10335 feet) do not affect the test results.
- [4] Hematocrit range: 10–70%.
- [5] Xylose may give falsely elevated results. DO NOT test blood glucose during or soon after xylose intake.
- [6] Tolazamide in therapeutic concentration may cause inaccurate results.
- [7] Acetaminophen and Paralidoxime Iodide (PAM) in concentrations higher than therapeutic concentration may cause inaccurate results.
- [8] Uric acid in concentrations higher than physiologic concentration may cause inaccurate results.

Performance Evaluation Data

Accuracy

A comparison study between the Wellion AMADEUS Meter and YSI 2300 Glucose Analyzer yields the following linear regression data:

Referring to EN ISO 15197:2015

Number of Readings: 600
Sample Range: 39 to 578 mg/dL (2,5 to 32,1 mmol/l)

Accuracy for glucose level <100 mg/dL (5,56 mmol/l)

Within ± 5 mg/dL (0,3 mmol/l)	Within ± 10 mg/dL (0,6 mmol/l)	Within ± 15 mg/dL (0,8 mmol/l)
120/174 (69,0%)	168/174 (96,6%)	174/174 (100%)

Accuracy for glucose level ≥100 mg/dL (5,56 mmol/l)

Within ± 5%	Within ± 10%	Within ± 15%
332/426 (60,5%)	419/426 (98,4%)	426/426 (100%)

Total within ± 15mg/dL (0,83 mmol/l) & ± 15% 600/600 (100%)
--

Precision

3 lots of Wellion AMADEUS Blood Glucose Test Strips were used for a within-run repeatability study. Glucose was added to venous blood in heparintubes in 5 different concentration. Blood glucose readings were recorded for 1 day resulting in 300 data points for each concentration. For evaluation results refer to the following table:

Repeatability Study
– According to EN ISO 15197:2015

Number of Readings:	300	300	300	300	300
Mean (mg/dL):	45	91	130	220	331
Mean (mmol/l):	2,5	5,1	7,2	12,2	18,4
S.D. (mg/dL)	2,4	2,6	3,5	6,3	6,8
S.D. (mmol/l)	0,1	0,1	0,2	0,4	0,4
C.V.%	NA	NA	2,7	2,8	2,0

3 lots of Wellion AMADEUS Test Strips were used for intermediate precision study. 3 levels of control material were prepared for glucose testing and 300 data points were obtained for each level, as shown in the following table:

Intermediate Precision Study
– According to EN ISO 15197:2015

Number of Readings:	300	300	300
Mean (mg/dL):	40	122	331
Mean (mmol/l)	2,2	6,8	18,4
S.D. (mg/dL)	3,1	3,2	5,0
S.D. (mmol/l)	0,2	0,2	0,3
C.V.%	NA	2,6	1,5

User Performance Study

A study evaluating glucose values from fingertip capillary blood samples obtained by 100 lay persons showed the following results: 100% within ± 15 mg/dL (0,83 mmol/L) of the YSI reference values at glucose concentration < 100 mg/dL (5,55 mmol/L) and 100% within ± 15% of the YSI reference values at glucose concentration ≥ 100 mg/dL (5,55 mmol/L).

Traceability

The system is calibrated using reference plasma values determined with a YSI analyzer. The YSI analyzer is calibrated using a series of standards traceable to NIST SRM917.

Serious incidents must be reported to the manufacturer and the national authorities.

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OBELIS S.A.